

**Anti-GCN4 tag antibody (HLENEVARLKK) [C11L34] (STJA0019301)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/SPR
Host / Source	Mouse
Reactivity	General

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.77 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**TARGET INFORMATION**

Immunogen Region	HLENEVARLKK
------------------	-------------

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--