

Anti-Kappa chain antibody [SAA0490] (STJA0019333)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/IF/IP/WB
Host / Source	Human
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.63 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-hFc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene Symbol	Igkc Uniprot ID= IGKC_MOUSE "
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--

