

Anti-LYZ antibody [SAA1157] (STJA0019666)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-hFc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	4069
Gene Symbol	LYZ
UniProt ID	LYSC_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--

