

Anti-Candida albicans ENO1/Enolase 1-strain SC5314/ATCC MYA-2876 & strain ATCC MYA-3404/T1 antibody [SAA1432] (STJA0019701)**GENERAL INFORMATION**

| | |
|---------------|-------------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC/SPR/WB |
| Host / Source | Human |
| Reactivity | C.albicans/C.tropicalis |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1.42 mg/mL |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG1 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|