

**Anti-MERS-CoV S2 Protein-isolate United Kingdom/H123990006/2012 antibody [G4]
(STJA0019729)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/IF/Neut/SPR
Host / Source	Human
Reactivity	MERS-CoV

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--