

Anti-Influenza A virus Neuraminidase/NA-A/environment/Shanghai/S1439/2013 & A/Shanghai/02/2013 & A/chicken/Jiangsu/S1460/2016 antibody [NA-63] (STJA0019738)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	Influenza A

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--