

**Anti-DENV-2 E-DII/Envelope Protein E antibody [2A10] (STJA0019774)****GENERAL INFORMATION**

|               |                      |
|---------------|----------------------|
| Product Type  | Primary antibodies   |
| Applications  | ELISA/Neut           |
| Host / Source | Mouse                |
| Reactivity    | DENV1-4/YFV/ZIKV/WNV |

**PRODUCT PROPERTIES**

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 4 mg/mL   |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Formulation         | 0.01M PBS, pH 7.4.  |
| Isotype             | IgG2a   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | <b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS</b> |
|------|---|