

## Anti-PLK1 antibody [3H10] (STJA0019789)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column
Formulation	0.01M PBS, pH 7.4.
Isotype	DARPin
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">5347</a>
Gene Symbol	<a href="#">PLK1</a>
UniProt ID	<a href="#">PLK1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--