

## Anti-CD19 scFv antibody [FMC63] (STJA0019809)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column
Formulation	0.01M PBS, pH 7.4.
Isotype	scFv-His
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">930</a>
Gene Symbol	<a href="#">CD19</a>
UniProt ID	<a href="#">CD19_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--