

Anti-SLC4A7 antibody [SAA1974] (STJA0020402)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.44 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Isotype	VHH-mFc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	9497
Gene Symbol	SLC4A7
UniProt ID	S4A7_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--