

Anti-Recombinant-FMC63 scFv antibody [SAA1996] (STJA0020527)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.43 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
UniProt ID	CD19_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--