

**Anti-CD56/NCAM1 antibody [SAA2007] (STJA0020554)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/WB           |
| Host / Source | Human              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Concentration       | 1 mg/mL                                                                                           |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Protein A/G purified from cell culture supernatant                                                |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                |
| Isotype             | IgG1                                                                                              |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">4684</a>        |
| Gene Symbol | <a href="#">NCAM1</a>       |
| UniProt ID  | <a href="#">NCAM1_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|