

## Anti-CD274/PD-L1/B7-H1 antibody [SAA2018] (STJA0020565)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-hFc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">29126</a>
Gene Symbol	<a href="#">CD274</a>
UniProt ID	<a href="#">PD1L1_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS

