

Anti-GAD67 antibody [L127/8] (STJA0021310)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC/WB |
| Host / Source | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG2a |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

ADDITIONAL INFORMATION

| | |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|