

Anti-Monkeypox virus/MPXV A29L antibody [SAA0286] (STJA0021600)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IB/ICC/IHC/Neut/WB
Host / Source	Human
Reactivity	Monkeypox virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.65 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--