

Anti-Mexican scorpion Beta-toxin Cl13/Cb3 antibody [11F] (STJA0021617)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut/SPR
Host / Source	Mouse
Reactivity	C.limpidus/C.baergi

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--