

**Anti-Staphylococcus aureus Alpha-toxin/hly/AT-strain NCTC 8325/PS 47 antibody [2A31]
(STJA0021628)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut/SPR
Host / Source	Mouse
Reactivity	S.aureus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--