

**Anti-GFP antibody [N86/20] (STJA0021658)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	A.victoria

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS</b>
------	---------------------------------------------------------------------------------------------------------------------------