

**Anti-DENV-2 Envelope protein E/EDIII domain-strain Thailand/16681/1984 antibody [mAb-61]
(STJA0021707)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IF/WB
Host / Source	Human
Reactivity	DENV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--