

Anti-DENV-1/2/3/4 Envelope protein E/EDIII A strand-strain Thailand/0348/1991 & Jamaica/1409/1983 & Brazil/97-11/1997 & Philippines/H87/1956 antibody [4E11#] (STJA0021712)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	DENV-4/DENV-2/DENV-1/DENV-3

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--