

Anti-DENV-1/3/4 Envelope protein E/EDIII domain-strain Nauru/West Pac/1974 & Philippines/H241/1956 antibody [2H12] (STJA0021713)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	DB/ELISA
Host / Source	Human
Reactivity	DENV-1/DENV-3/DENV-4

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--