

Anti-Suid herpesvirus 1/SuHV-1 UL42/PAP-strain Kaplan antibody [2G11] (STJA0021732)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IF/IP/WB
Host / Source	Mouse
Reactivity	SuHV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--