

Anti-IL21 antibody [SAA2070] (STJA0021757)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/WB
	Host / Source	Human
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	1 mg/mL
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant
	Formulation	0.01M PBS, pH 7.4.
	Isotype	IgG1
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months
TARGET INFORMATION		
	Gene ID	59067
	Gene Symbol	IL21
	UniProt ID	IL21_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS