

Anti-TREM2 antibody [aI002c] (STJA0021761)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Agonist/DB/EMSA/FISH/ICC/IF/IP/SPR
	Host / Source	Mouse
	Reactivity	Human/Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	5.68 mg/mL
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant
	Isotype	IgG1
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months
TARGET INFORMATION		
	Gene ID	83433
	Gene Symbol	Trem2
	UniProt ID	TREM2_MOUSE
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS