

Anti-Norovirus/Norwalk virus Major capsid protein/VP1-Hu/GII4/Farmington Hills/2004/USA antibody [A1431] (STJA0021767)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Human
Reactivity	Norovirus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.13 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--