

Anti-B2M/Beta 2-MG antibody [1A053] (STJA0021773)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/GICA/CLIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	PBS buffer, pH 7.4, containing 0.1% NaN3 as preservative.
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	567
Gene Symbol	B2M
UniProt ID	B2MG_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--