

## Anti-ZBTB7B antibody [C206] (STJA0021779)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IF/IP        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Concentration       | 1 mg/mL                                                                                           |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Protein A/G purified from cell culture supernatant                                                |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                |
| Isotype             | IgG2a                                                                                             |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">51043</a>       |
| Gene Symbol | <a href="#">ZBTB7B</a>      |
| UniProt ID  | <a href="#">ZBT7B_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|