

Anti-CD5/Ly-1 antibody [SAA2083] (STJA0021869)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | Human |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | VHH-hFc |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 12507 |
| Gene Symbol | Cd5 |
| UniProt ID | CD5_MOUSE |

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS**

