

Anti-BTNL2/BTL-II antibody (Asp27-Gly453) [1A066] (STJA0021904)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	547431
Gene Symbol	Btl2
UniProt ID	BTNL2_MOUSE
Immunogen Region	Asp27-Gly453

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--