

**Anti-ENO1/Alpha-enolase antibody (Met1-Lys434) [1A067] (STJA0021905)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**TARGET INFORMATION**

Gene ID	<a href="#">2023</a>
Gene Symbol	<a href="#">ENO1</a>
UniProt ID	<a href="#">ENOA_HUMAN</a>
Immunogen Region	Met1-Lys434

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--