

Anti-Agrobacterium sp aroA/EPSPS antibody (His3-Ala455) [1A079] (STJA0021917)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Agrobacterium sp

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Immunogen Region	His3-Ala455
------------------	-------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--