

## Anti-KLB/Beta-klotho antibody (Lys2-Pro996) [1A138] (STJA0021965)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.3 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">152831</a>
Gene Symbol	<a href="#">KLB</a>
UniProt ID	<a href="#">KLOTB_HUMAN</a>
Immunogen Region	Lys2-Pro996

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--