

## Anti-CD38 antibody (Val43-Ile300) [1A166] (STJA0022009)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">952</a>
Gene Symbol	<a href="#">CD38</a>
UniProt ID	<a href="#">CD38_HUMAN</a>
Immunogen Region	Val43-Ile300

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--