

Anti-Bacteriophage M13 Capsid protein G8P antibody [1A191] (STJA0022024)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Phage M13

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--