

## Anti-GST-tag antibody [1A214] (STJA0022048)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/WB           |
| Host / Source | Mouse              |
| Reactivity    | S.japonicum        |

### PRODUCT PROPERTIES

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Concentration       | 1mg/mL                                                                                            |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Protein A/G purified from cell culture supernatant                                                |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                |
| Isotype             | IgG1                                                                                              |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

### ADDITIONAL INFORMATION

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|