

Anti-Clostridium Difficile Toxin B/tcdB antibody [1A252] (STJA0022084)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	GICA/IFA
Host / Source	Mouse
Reactivity	C.difficile

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--