

Anti-ALB/Albumin/HSA antibody [1A390] (STJA0022222)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/TIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.62 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.02M PBS, pH 7.4, containing 0.1% Preservatives
Isotype	IgG
Storage Instruction	Store at +4°C for twelve months from the date of receipt

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	---