

Anti-IgG Kappa Chain antibody [SAA0278] (STJA0022297)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	VHH-mFc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene Symbol	IGKC Uniprot ID= IGKC_HUMAN "
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--