

## Anti-GFP Polyclonal antibody (Met1-Lys238) (STJA0022380)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Rabbit
Reactivity	GFP-Tag

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2.46 mg/mL
Conjugation	Unconjugated
Purification	Purified by antigen affinity column
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Immunogen Region	Met1-Lys238
------------------	-------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--