

Anti-6*His Tag antibody [1A056] {HRP} (STJA0022404)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	General

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.48 mg/mL
Conjugation	HRP
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--