

Anti-FOXP3 antibody (Ser304-Pro431) [1A208] {APC} (STJA0022424)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	APC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	50943
Gene Symbol	FOXP3
UniProt ID	FOXP3_HUMAN
Immunogen Region	Ser304-Pro431

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--