

## Anti-FOXP3 antibody (Ser304-Pro431) [1A209] {APC} (STJA0022425)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC/WB        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | APC   |
| Purification        | Affinity chromatography                                     |
| Formulation         | 0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.         |
| Isotype             | IgG1  |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">50943</a>       |
| Gene Symbol      | <a href="#">FOXP3</a>       |
| UniProt ID       | <a href="#">FOXP3_HUMAN</a> |
| Immunogen Region | Ser304-Pro431               |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|