

## Anti-FOXP3 antibody (Ser304-Pro431) [1A211] {APC} (STJA0022427)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	APC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

### TARGET INFORMATION

Gene ID	<a href="#">50943</a>
Gene Symbol	<a href="#">FOXP3</a>
UniProt ID	<a href="#">FOXP3_HUMAN</a>
Immunogen Region	Ser304-Pro431

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--