

Anti-FOXP3 antibody (Ser304-Pro431) [1A209] {PE} (STJA0022429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Purification	Affinity chromatography
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

TARGET INFORMATION

Gene ID	50943
Gene Symbol	FOXP3
UniProt ID	FOXP3_HUMAN
Immunogen Region	Ser304-Pro431

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------