

Anti-6*His Tag antibody [1A056] {APC} (STJA0022434)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Mouse
Reactivity	General

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	APC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--