

Anti-EGFP antibody [1A190] {HRP} (STJA0022439)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/WB
Host / Source	Mouse
Reactivity	Water jellyfish

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--