

**Anti-Bacteriophage M13 Capsid protein G8P antibody [1A191] {FITC} (STJA0022440)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA              |
| Host / Source | Mouse              |
| Reactivity    | Bacteriophage M13  |

**PRODUCT PROPERTIES**

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | FITC  |
| Purification        | Affinity chromatography                                     |
| Formulation         | 0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.         |
| Isotype             | IgG1  |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

**ADDITIONAL INFORMATION**

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|