

Human IgG1 Isotype Control [REA293] {Unconjugated} (STJA0022494)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	3.15 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS**

