

Mouse IgM Isotype Control [11E10#] {Unconjugated} (STJA0022541)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	3.4 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	30 mM Histidine, pH5.8, 10% sucrose, 0.02% Tween80
Isotype	IgM
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS

