

Monkey IgG Isotype Control [NA] {Unconjugated} (STJA0022626)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	4.2 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--